



Introduction to Beam Diagnostics

Ulrich Raich
CERN BE - BI
(Beam Instrumentation)



A few depicted examples




- Introduction
- Beam presence: Scintillating screens (LHC)
- Intensity measurement, Faraday Cup and Transformer (Linac-4)
- Transverse Profile measurement, wire scanner & wire grids (PSB & PS)
- Emittance measurement
 - Slit and Grid (Linac-4, 3 MeV line)
 - Emittance measurement line (Linac-2)
 - Longitudinal Emittance measurement (Linac-2)
- Trajectory measurement (LHC and PS) using Beam Position Monitors (BPMs)
- Longitudinal phase space: Tomoscope (PS) using a wall current monitor
- Tune measurement (SPS) using BPMs
- Losses: Beam Loss Monitors (BLMs) (LHC)



Introduction



An accelerator can never be better than the instruments measuring its





Diagnostic devices and quantity measured




| Instrument | Physical Effect | Measured Quantity | Effect on beam |
|----------------------------|---------------------------------------|-----------------------------------|--|
| Faraday Cup | Charge collection | Intensity | Destructive |
| Current Transformer | Magnetic field | Intensity | Non destructive |
| Wall current monitor | Image Current | Intensity Longitudinal beam shape | Non destructive |
| Pick-up | Electric/magnetic field | Position, Tune | Non destructive |
| Secondary emission monitor | Secondary electron emission | Transverse size/shape, emittance | Disturbing, can be destructive at low energies |
| Wire Scanner | Secondary particle creation | Transverse size/shape | Slightly disturbing |
| Scintillator screen | Atomic excitation with light emission | Transverse size/shape (position) | Destructive |
| Residual Gas monitor | Ionization | Transverse size/shape | Non destructive |



A beam parameter measurement





U. Raich, CERN School of Accelerators, Chavannes 2013/14



Beam Presence Scintillating Screens



Method already applied in cosmic ray experiments


- Very simple
- Very convincing

Needed:

- Scintillating Material
- TV camera
- In/out mechanism

Problems:

- Radiation resistance of TV camera
- Heating of screen (absorption of beam energy)
- Evacuation of electric charges






Screen mechanism



Screen with graticule









Results from TV Frame grabber






First full turn as seen by the BTV 10/9/2008



Un-captured beam sweeps through he dump line




 For further evaluation the video signal is digitized, read-out and treated by program



Beam Intensity Layout of a Faraday Cup



- Electrode: 1 mm stainless steel
- Only low energy particles can be measured
- Very low intensities (down to 1 pA) can be measured
- Creation of secondary electrons of low energy (below 20 eV)
- Repelling electrode with some
 100 V polarisation voltage pushes
 secondary electrons back onto the electrode






Faraday Cup










Electro-static Field in Faraday Cup



In order to keep secondary electrons with the cup a repelling voltage is applied to the polarization electrode

Since the electrons have energies of less than 20 eV some 100V repelling voltage is sufficient






Energy of secondary emission electrons

I(µA)



- With increasing repelling voltage the electrons do not escape the Faraday Cup any more and the current measured stays stable.
- At 40V and above no decrease in the Cup current is observed any more





Intensity Principle of a fast current transformer






Diagram by H. Jakob




Magnetic shielding



Transformer steel (µ2) Soft iron (µ1)

Permalloy (µ3)






- Shield should extend along the vacuum chamber length > diameter of opening
- Shield should be symmetrical to the beam axis
- Air gaps must be avoided especially along the beam axis
- Shield should have highest µ possible but should not saturate



Calibration of AC current transformers






- The transformer is calibrated with a very precise current source
- The calibration signal is injected into a separate calibration winding
- A calibration procedure executed before the running period
- A calibration pulse after the beam pulse measured with the beam signal



Profile measurements



Secondary emission grids (SEMgrids)



The ejected electrons are taken away by polarization voltage

When the beam passes secondary electrons are ejected from the wires


The current flowing back onto the wires is measured





Profiles from SEMgrids





Projection of charge density projected to x or y axis is Measured

One amplifier/ADC per wire Large dynamic range


Resolution is given by wire distance

Used only in transfer lines


Waiting for new measurement timing !!!



Wire Scanners



A thin wire is quickly moved across the beam Secondary particle shower is detected outside the vacuum chamber on a scintillator/photo-multiplier assembly Position and photo-multiplier signal are recorded simultaneously








Wire scanner profile





High speed needed because of heating.

Adiabatic damping


Current increase due to speed increase

Speeds of up to 20m/s => 200g acceleration



Emittance measurements





Design by E. Bravin

A beam is made of many, many particles, each one of these particles is moving with a given velocity. Most of the velocity vector of a single particle is parallel to the direction of the beam as a whole (s). There is however a smaller component of the particles velocity which is perpendicular to it (x or y).


$$\vec{v}_{particle} = v_s \hat{u}_s + v_x \hat{u}_x + v_y \hat{u}_y$$



Emittance measurements




- If for each beam particle we plot its position and its transverse angle we get a particle distribution whose boundary is an usually ellipse.
- The projection onto the x axis is the beam size





The slit and grid method

- If we place a slit into the beam we cut out a small vertical slice of phase space
- Converting the angles into position through a drift space allows to reconstruct the angular distribution at the position defined by the slit







Transforming angular distribution to profile




Influence of a drift space

- When moving through a drift space the angles don't change (horizontal move in phase space)
- When moving through a quadrupole the position does not change but the angle does (vertical move in phase space)






Influence of a quadrupole






The Slit Method










Emittance Meter








Transverse Emittance Measurement





Slit and grid phase space scanner

L-shaped 0.1mm slit moves under 45 degrees

Slit and grids move independently Positioning precision: 50 µm Movement PLC controlled

Slit and grids mounted in 2 independent vacuum boxes which can be separated


Horizontal and vertical SEMGrid

- wire distance .75 mm
- 30 signal wires
- readout with home built 36 channel 250 kHz ADC
- time resolved profiles



Emittance plot Solenoid








Single pulse emittance measurement






Quadrupole

Kickers



slit



Quadrupole


SEMgrid





Transformation in Phase Space








Result of single pulse emittance measurement






Waiting for new acquisition...



Longitudinal Emittance measurement






- second slit selects energy slice
- first kicker sweep phase space over all energies
- buncher rotates energy slice in phase space
- at second spectrometer the phase distribution is transformed into an energy distribution analyzed by the second spectrometer
- second kicker corrects for first kick

Ε


at slit





Trajectory and Orbit Measurements








The PS, a universal machine




All beams pass through the PS
Different particle types
Different beam characteristics
Concept of a super cycle





The Supercycle







Position Measurements



BPM signals sampled at 120 MHz


Red: The sum signal

Green: The difference signal

Procedure:

Produce integration gates and baseline signals
Baseline correct both signals
Integrate sum and difference signals and store results in memory
Take external timing events into account e.g. harmonic number change, γ-transition etc.








Baseline restoration






Low pass filter the signal to get an estimate of the base line Add this to the original signal




Trajectory measurements in circular machines



Needs integration gate
Can be rather tricky
Distance between bunches
changes with acceleration
Number of bunches
may change






U. Raich, CERN School of Accelerators, Chavannes 2013/14



Beams in the PS










RF Gymnastics










- Bunch splitting or recombination
 One RF frequency is gradually
 decreased while the other one
 is increased
- The gate generator must be synchronized



Trajectory readout electronics








Tune measurements




- When the beam is displaced (e.g. at injection or with a deliberate kick, it starts to oscillate around its nominal orbit (betatron oscillations)
- Measure the trajectory
- Fit a sine curve to it
- Follow it during one revolution





The kicker


The Sensors




Shoebox pick-up with linear cut













Kicker + 1 pick-up



- Measures only non-integral part of Q
- Measure a beam position at each revolution
- Fourier transform BPM signal
- Search peak in Fourier Spectrum






Fourier transform of pick-up signal



Q-Measurement Results








Direct Diode Detection Base Band Q measurement





Diode Detectors convert spikes to saw-tooth waveform

Signal is connected to differential amplifier to cut out DC level

Filter eliminates most of the revolution frequency content


Output amplifier brings the signal level to amplitudes suitable for long distance transmission




BBQ Results from CERN SPS



Results from Sampling










Tune feedback at the LHC







Computed Tomography (CT)



Principle of Tomography:


- Take many 2-dimensional Images at different angles
- Reconstruct a 3-dimensional picture using mathematical techniques (Algebraic Reconstruction Technique, ART)








The reconstruction











Produce many projections of the object to be reconstructed


Back project and overlay the "projection rays" Project the backprojected object and calculate the difference


Iteratively backproject the differences to reconstruct the original object



Some CT resuluts










Computed Tomography and Accelerators




RF voltage

Restoring force for nonsynchronous particle

Longitudinal phase space


Projection onto Φ axis corresponds to bunch profile





Reconstructed Longitudinal Phase Space







Courtesy S. Hancock



Bunch Splitting










Stored Beam Energies





| Quench Levels | Units | Tevatron | RHIC | HERA | LHC |
|-----------------------------|----------------------|-----------------------|-----------------------|---|-----------------------|
| Instant loss (0.01 - 10 ms) | [J/cm ³] | 4.5 10 ⁻⁰³ | 1.8 10 ⁻⁰² | 2.1 10 ⁻⁰³ - 6.6 10 ⁻⁰³ | 8.7 10 ⁻⁰⁴ |
| Steady loss (> 100 s) | [W/cm ³] | 7.5 10 ⁻⁰² | 7.5 10 ⁻⁰² | | 5.3 10 ⁻⁰³ |



Beam power in the LHC





The Linac beam (160 mA, 200µs, 50 MeV, 1Hz) is enough to burn a hole into the vacuum chamber

What about the LHC beam: 2808 bunches of 15*10¹¹ particles at 7 TeV?

1 bunch corresponds to a 5 kg bullet at 800 km/h




Beam Loss Monitor Types



- Design criteria: Signal speed and robustness
- Dynamic range (> 10⁹) limited by leakage current through insulator ceramics (lower) and saturation due to space charge (upper limit).

Secondary Emission Monitor (SEM):

- Length 10 cm
- $P < 10^{-7} bar$
- ~ 30000 times smaller gain



Ionization chamber:


- N₂ gas filling at 100 mbar over-pressure
- Length 50 cm
- Sensitive volume 1.5 l
- lon collection time 85 μs
- Both monitors:
 - Parallel electrodes (Al, SEM:
 Ti) separated by 0.5 cm
 - Low pass filter at the HV input





Quench levels








Industrial production of chambers





Beam loss must be measured all around the ring => 4000 sensors!



Conclusions



- Beam diagnostics is a very wide field where many different competences are needed
 - Machine physics
 - Electronics
 - Computing
 - Mechanics
- The instruments are the eyes with which we observe the beam
- The beam can never be adjusted with higher precision than what can be measured